

EVA LR-BV PLUS

FEATURES AND BENEFITS

- Temperature range +12°C to +60°C (according to models).
- Compatible with all PLUS range viscometer.
- Rapid temperature setting of 600 ml and 150 ml beakers according to ISO 2555.



ACTIVITY DOMAINS

-  Agro-Alimentaire
-  Cosmetics pharmaceuticals
-  Paint / ink / coatings
-  Chemical / petroleum products
-  Car industry
-  Teaching

TEMPERATURE CONTROL DEVICES AVAILABLE

Compatible only with bayonet coupling device

Designation	Part Number	Temperature range	Type	Compatible device	Compatible Measuring system	Size (pxlxh) in mm	Weight
EVA LR PLUS ^{a)}	N950006	+15°C to +45°C	Peltier	All	ASTM (L and R)	610x340x650	15Kg
EVA BV PLUS ^{b)}	N950300	+12°C to +60°C	Peltier		MS-BV		

a) Usable with 600 ml beaker.

b) Usable with 150 ml beaker. Not compatible with all instruments version LR