

EVA MS-R PLUS

FEATURES AND BENEFITS

- Temperature range +12°C to +60°C (according to models).
- Compatible from RM 100 PLUS.
- Can accommodate two instruments.
- Support for 9 simultaneous samples.



ACTIVITY DOMAINS



Food Industry



Cosmetics pharmaceuticals



Paint / ink / coatings



Chemical / petroleum products



Car industry



Building materials

TEMPERATURE CONTROL DEVICES AVAILABLE

(Not compatible with B-ONE and LR version)

Compatible only with bayonet coupling device.

Designation	Part Number	Temperature range	Type	Compatible device	Compatible Measuring system	Size (pxlxh) in mm	Weight
EVA MS-R PLUS 9 Holes	N950200	+17°C to +45°C	Peltier	RM 100 ^{a)} RM 200 DSR 500	MS-R	610x340x653	15Kg
EVA MS-R PLUS 9 Holes	N950210	+12°C to +60°C	Peltier				

^{a)} Usable with RM100 Portable and RM100 i